




**CONEFIT**  
Implant System






# Table of contents




## INTRODUCTION

About the company .....	3
Technology and quality .....	4
Implant surface finish .....	6
Applied raw materials .....	7
Packaging .....	8



## CONEFIT IMPLANT SYSTEM

The distinctive characteristics of the Conefit system .....	14
The applicational fields of Conefit .....	16
Bi-phase Conefit implants .....	18
Single-phase, two-piece Conefit implants .....	20
Single-phase, one-piece „A“ type Conefit implants .....	22
Single-phase, one-piece „B“ type Conefit implants .....	24
Drilling protocol .....	26




## CONEFIT ABUTMENT SYSTEM

Conefit abutments .....	28
-------------------------	----



## INSTUMENTS

Conefit instrument kits .....	40
Surgical drills .....	44
Ratchet torque wrench .....	46





# About the company

**BIONIKA Medline** Orvostechikai Kft. Is a member of the Swedish-Hungarian group of companies. It was founded in 1989. The owners of the company are Swedish and Hungarian citizens.

We have more than 30-year-experience in the field of medical instruments and implant development, production and trade. BIONIKA as a researcher, developer, manufacturer and distributor is present in dentistry, oral surgery, traumatology, orthopedics and rehabilitation in the medical-professional areas.

According to our objective and perception, we attach great importance to the word „BIONIKA”, which marks a scientific thinking on the boundaries of biology, technology and electronics that combines these three areas in our researching and developing work.

**Clinical and technological experiences:** The continuous process, combination and utilization of clinical and technological experiences in development contributes to our success, up to the production base. Here you will find the best solutions and constructions suited to customer needs, which are under continuous development.


**Development:** The owners of BIONIKA put great emphasis on continuous product and technological research and development. Our products are developed in close collaboration with doctors and engineers, enabling us to ensure the world-class quality and practical utilization.

**Quality:** The quality of the products expected by our customers is guaranteed by design, manufacturing and quality management according to the harmonized European Union laws. The BIONIKA Medline Kft. is operated according to the EN ISO 9001 and the EN ISO 13485 quality management system. Our products are provided with CE marks.

**Guarantee:** After inserting the implant - the risk of the ossification process is assumed by BIONIKA, independently of cause and effect relationship – exchange guarantee is ensured within one year after the purchase. Otherwise, we provide a long-term, 10-year guarantee for our products.



**DENTISTRY**



**ORAL SURGERY**



**TRAUMATOLOGY**



**ORTHOPEDICS**

# Technology

BIONIKA Medline Kft. has more than 30 years of experience in the development and production of dental implants, dental insertion instruments and stomatological parts. During this time more than 40 types of implant systems have been developed and are being manufactured to date, including insertion instruments. Some of these parts have been developed for their own marketing in accordance with their own market needs.

Other systems – in cooperation with independent medical groups – are made to order, mainly developed and manufactured for foreign markets. (These are sold by the customers under their own brand name).

Our partners can choose from approximately 20.000 different parts of different sizes and shapes. Our manufacturing technology is flexible, we can quickly move from one component to another, and we are able to fulfill thousands of orders with a short turnaround time. This area requires high precision production (in some cases it is necessary to hold 2-5 mm tolerances).

All the technological operations we carry out are from manufacturing, surface design, packaging. Our products are CE marked and the production process is under strict quality management system. Biocompatible materials are the most important raw materials for dental, oral surgery, traumatology and orthopedic medical implants. Because relatively small series of customized solutions are required, they require fast programmable CNC machining technology.

Accordingly, we have molded CNC machining centers and Swiss type longitudinal machining centers. For machining more complex surfaces, an industrial 5-axis CNC center is used with CAD-CAM system support. Our machines are equipped not only with fixed, but also with propelled cutting instrument units, with which we can perform more complex spatial geometrical machining. As a complementary technology, we have sandblasting, polishing titanium coloring and sterilization equipments. The production of custom prosthetic components for dental applications is supported by the BIONIKA Milling center.

## Our Partners



CITIZEN

SIEMENS

SAASCO  
medical devices Med

phv  
PHARMAVALID Kft.

FRESH-  
DENTAL

OGYÉI  
Országos Gyógyszerészeti  
és Élelmezés-egészségügyi Intézet



formlabs

bay

graphIT

NCT  
CONTROL DRIVES MOTORS  
KEEP MOVING

EIKI  
MSZ EN ISO 9001:2009  
MSZ EN ISO 13485:2004

AGROSTER  
Besugárzó Rt.


MEDICOR



MISKOLCI  
EGYETEM  
UNIVERSITY OF MISKOLC

DrJuice Co.  
A teremtő megalkotta, mi tökéletesítettük

mta ttk



# Quality management and guarantee

The quality of the products is guaranteed by design, manufacturing and quality management according to the harmonized European Union laws. The BIONIKA Medline Kft. is operated according to the EN ISO 9001 and the EN ISO 13485 quality management system. Our products are provided with CE marks, which was ensured by EMKI and QT-CERT.

We provide a long-term, 10-year guarantee for our products. After inserting the implant, reducing the medical risk of the ossification process, independently of cause and effect relationship – exchange guarantee is ensured within one year after the purchase for the dropped and fallen out implants.



BIONIKA Medline Kft. has always paid close attention to quality and reliability during its more than 30 years of existence. The Bisnode certificate is proof of our reliability and stability. BIONIKA received „Triple A” Bisnode qualification from 2016 to 2022.

At the Hungarian market only 0.63% of the companies have the AAA Bisnode rating. The financial risk of establishing business relationship with these business associations is extremely low - source: dnb.hu


# Superclean implant surface

**BIONIKA demonstrates the best qualities of Grade 4 titanium used in implant manufacturing for dental implantology according to the ISO 5832-2 ASTM F67 standard. Due to its adequate purity the biocompatibility is exceptionally good as well as it is provided with exceptional solidity.**

Initially, we and other implant manufacturers preferred the higher purity titanium but due to solidity reasons nowadays almost every implant is made of Grade 4 or other alloyed titanium in the world.

In all cases of implant abutments, alloyed, high strength Grade 5 titanium is applied according to the ISO 5832-4 ASTM F136 standard. The titanium applied according to the standard is provided with exceptional biocompatibility, it is almost risk-free.

Almost all professionals see that the implantation success is best determined by the implantologist's practice, as well as surgical conditions, carefully maintained hygiene, and patient abilities.




Energy dispersive X-ray spectrometric elemental analysis of Bionika implants\*

**The main steps of our BioTiS surface finish technology:**

- Chemical, mechanical surface cleaning and surface dewing
- Special ultrasonic cleaning, surface cleaning and sterilization
- Transformation of surface structure by acidification process
- Multi-stage dehumidification, cleaning
- Electrochemical surface modification
- sterilization
- Surface finish in physiological solution

These technological steps are always carried out under sterile conditions.

The final packaging of the implants is four-layered. The packaging is carried out in a sterile cabin. Final sterility is assured by an accredited 20 Rad gamma sterilization procedure.



\*Source: FOGSurgical SZEMLE, year 106. No. 4.. 2013. 135-143

Bionika implant electron microscope i\*



# Applied raw materials



## Titanium grade 4

### Chemical composition

Elements	Threshold limit of constituents(%)
O	0,4 max.
Fe	0,3 max.
C	0,1 max.
N	0,05 max.
H	0,0125 max.
Ti	>99% / balance

### Mechanical properties

solidity	680 MPa min.
dilation	10 %

According to the **ISO 5832-2** standard.

## Titanium Grade 5

### Chemical composition

Elements	Threshold limit of constituents(%)
Al	5,5-6,75 max.
V	3,5-4,5 max.
Fe	0,3 max.
O	0,2 max.
C	0,08 max.
N	0,05 max.
H	0,015 max.
Ti	balance

### Mechanical properties

solidity	860 MPa min.
dilation	10 %

According to the **ISO 5832-3** standard.

## CoCr

### Chemical composition

Elements	Threshold limit of constituents(%)
C	0,1 max.
Si	1,0 max.
Mn	1,0 max.
P	0,005 max.
S	0,005 max.
Cr	30, 0 max.
Mo	7,0 max.
Ni	1,0 max.
Co	-
N	0,2250 max.

### Mechanical properties

solidity	1240,00 MPa min.
elongation limit	900,00 min.
elongation at break	18,00 min.
fracture contraction	23,00 min.

According to the **ISO 5832-4** standard.

## Plastics

**POM** (polyoxymethylene) : Thermoplastic synthetic plastic, Excellent properties eg: high hardness, low wear, good flexibility, little absorbing ability. Density: 1.41 g / cm<sup>3</sup>. elongation at break: min. 30% Current Voltage: min. 65 Mpa. Its color is white.

**PEEK** (polyether ether ketone) :High heat-resistant plastic, suitable for all conventional sterilization methods (steam, dry heat, ethylene oxide, gamma radiation). Density: 1.30 1.41 g / cm<sup>3</sup> Tensile strength: 115 Mpa. elongation at break: min. 17% Its colour is natural brownish gray.

# Conefit Packaging




10 - piece collection box



vial in sterile foil



removing the foil




sterile vial




removing the locking plug  
(1)




removing the locking plug  
(2)



using the implant key  
driver



removing the implant



## Collection box

Depending on the order quantities, collection boxes with 5 and 10 pieces are applied.

# Conefit Packaging




According to our endeavor, the raw materials used in packaging are almost 100% naturally occurring, environmentally-friendly, natural materials.

aluminum · titanium · glass · caoutchouc · paper

## Vial

The first layer of the packaging is a transparent vial, which ensures the total, null colony-forming unit sterility. The locking plug of the vial holds the implant, the multifunctional implant head and the locking screw.




## Paper box

The outer layer of the packaging is a paper box with a high density, which is for the physical safety. Every paper box is provided with colour-coded labels according to the different platform- diameters. The colour of the packaging is adjusted to this method.








## The sectional image of the **Conefit** packaging and its accessories

The inner layer of the packaging is the implant holding poor box. The poor box is held by the vial locking plug and with this they can be removed from the vial. The implant locking screw can be found in the vial locking plug as well.




# Conefit product labels and their notation

Differential platform diameters by colour and diameter:


-  3,5 mm - grey
-  3,75 mm - yellow
-  4,0 mm - purple
-  4,5 mm - blue
-  5,0 mm - green

Information supplied by the three product labels to the outer packaging of the Implant System:


The side of the box:

















Top of the box:




The back of the box:



## Explanation of label codes:

<ul style="list-style-type: none"> <li> Diameter</li> <li> Implant length</li> <li> Item number</li> <li> Serial number</li> <li> Production time</li> </ul>	<ul style="list-style-type: none"> <li> Expiration date</li> <li> It is forbidden to use in the case of damaged packaging.</li> <li> Reuse is forbidden!</li> <li> Sterilized with gamma rays</li> <li> sterilized with steam or dry heat</li> </ul>	<ul style="list-style-type: none"> <li> Non-sterile product in the package</li> <li> Read the usage guide!</li> <li> Certification company code</li> <li> Manufacturer's name and Contact</li> </ul>
---	---	--


**CONEFIT**  
Implants





## IMPLANT SYSTEM | TABLE OF CONTENTS



Distinctive characteristics of the Conefit implant system ...	14
The applicational fields of Conefit.....	16
Sizes available of the implants .....	18
Drilling protocol .....	26



# Distinctive characteristics of the **Conefit** implant system

The creation and development of the CONEFIT dental implant system was started by the BIONIKA Engineer Office in 1992 with the medical support of the Central Institute of Stomatology. The implantation experience and research of the Oral Surgery Department have led

to the recognition that the application of a self-locking, conical, high pitch titanium screw is a biomechanically advantageous solution for dental implantation. It is possible to have single and bi- phase solutions from surgical point of view.



## **Cortical level**

The chances of the implant persistency are significantly improved by inserting the upper edge of the implant at the cortical level or below (from 0,5 to 1 mm) which is proven by surgical experiences on the long run.

## **Platform switching**

The diameter of the abutment is smaller than the outer part of the implant which is connected to the bone. This way the soft tissue closes upon the ingoing aperture and the implant's connecting surface of the bone nest, covering, insulating that like as a sealing ring.

## **Spirally microstriated cortical surface**

The multi- paragraphed microstriated spiral surface can function as a significant weight bearing element connected to the cortical and ensuring micromotor-free condition for the fast inserting. This self-closing thread structure due to the cycloid cord thread geometry helps facilitating dynamic force convection and ensures micromotor-free condition for fast insertion.

## **Anatomical root form**

The Conefit follows the form of the anatomical tooth root. Due to the conicity, high thread pitch, high thread deepness, self-closing and self-tapping shaping of the implant screw thread, it has a bone-compacting effect. Thus, there is an option for immediate loadness of the implant as needed as a result of the high primer stability.

## **Rounded implant end**

It helps facilitating the minor direction changes when inserting the implant.



# Types of the **Conefit** system

## Bi-phase implant

The usage of the bi-phase Conefit system allows a wide range of variability. Our bi-phase implants are available in more than 40 kinds of geometry and in large range of sizes of abutments for plastic, zirconium, titanium and cobalt chromium-based toothworks.



**Connection**  
The head is provided with a **hexagonal** key aperture and a 60 degree centering morse cone, which is stably fit into the joining part of the abutment.


## Single-phase , one-piece implant

The Conefit single-phase, two-piece implant is manufactured in a collar design. It has basic features, besides that the main feature of it is the neck length adjusted to the average thickness of the mucous membrane.




## Single-phase , two-piece implant

The Conefit single-phase, one-piece Implant with anatomical head for fixing toothworks. "A" and "B" type.



# The applicational fields of the **Conefit** Implant system




## In the case of one tooth deficiency

In this case of the replacement of a tooth, we do not have to grind two healthy teeth for bridge replacement, but inserting an implant, then we need to glue a crown in the same way as the traditional one.

## In the case of end of line tooth deficit(s)

In this case, in the absence of a pillar tooth, we are not able to make a fixed replacement (bridge). With the implantation of at least two implants, you are already make the (fixed) bridge replacement.



## Removeable denture


## In the case of total tooth deficit

In this case the patient has no tooth, complete tooth augmentation can be performed. In this case, there one solution is the removable denture: 2-4 implants are implanted, these will be the fixation for the removable tooth.

This brings a tremendous quality of life to the patient, as this way the denture will be very stable, so it can be used in chewing and speaking outright. There are two solutions possible in this case: ball head or locator head abutments can be applied.

With the implantation of several 6-8 implants, it is possible to make full fixation (round bridge) augmentation, which is both functional and aesthetically close to the natural teeth.

## Screw-retained fixed dental prosthetics



**Optimum Concept**

## Optimum Concept

**All-on-4® type - Economical Solution**

**The Optimum Concept provides great stability, with only four implants being implanted.**


- The temporary denture can be inserted on the day of surgery.
- Immediate improvement in function, speech and aesthetically.
- Treatment times are shorter and costs can be lower than conventional implant placement modes.
- Tilt rear implants can be fixed better into the front bone. This promotes prosthesis support.

## Safe Concept

**All-on-6® type - For extra stability**

**The stability of the toothwork can be increased with the Safe Concept. It is exceptionally advantageous in the case of extra chewing ability.**


- The usage of oblique head implants allows longer implants to be used, avoiding the nerve canal.
- The usage of longer implants allows the bone and the implant to touch on a larger surface, thus making bone augmentation avoidable.
- Favorable bone level for tilted and axial implants.
- High remaining chances.



**Safe Concept**


# Conefit implant system

## Sizes available of the Conefit bi-phase implants








The Conefit bi-phase implant system consists of implants with five different diameters. Each diameter implant is available in five lengths to be able to find the right solution for each situation. The  $\varnothing 3.5$  diameter implant is exceptionally suitable for thinner bone than average for holding the toothworks in the long run. The diameter of  $\varnothing 3.75$ ,  $\varnothing 4.0$  and  $\varnothing 4.5$  is beneficial for the average bone structure, while the usage of  $\varnothing 5.0$  implants is advantageous for larger bone supply than average.

All our Conefit implants are made of a homogeneous, high-strength titanium alloy.








**Ø 3.5**

inserting length (L):

				
8 mm	10 mm	12 mm	14 mm	16 mm

**Ø 3.75**

inserting length (L):

				
8 mm	10 mm	12 mm	14 mm	16 mm

**Ø 4.0**

inserting length (L):

				
8 mm	10 mm	12 mm	14 mm	16 mm

**Ø 4.5**

inserting length (L):

				
8 mm	10 mm	12 mm	14 mm	16 mm


**Ø 5.0**

inserting length (L):

				
8 mm	10 mm	12 mm	14 mm	16 mm


# Conefit implant system

## Sizes available of the Conefit single-phase, two-piece implants



The Conefit single-phase, two piece implant system consists of implants with five different diameters. Each diameter implant is available in five lengths to be able to find the right solution for each situation. The Ø3.5 diameter implant is exceptionally suitable for thinner bone than average for holding the toothworks in the long run. The diameter of Ø3.75, Ø4.0 and Ø4.5 is beneficial for the average bone structure, while the usage of Ø5.0 implants is advantageous for larger bone supply than average.

All our Conefit implants are made of a homogeneous, high-strength titanium alloy.








**Ø 3.5**

inserting length (L):

				
8 mm	10 mm	12 mm	14 mm	16 mm

**Ø 3.75**

inserting length (L):

				
8 mm	10 mm	12 mm	14 mm	16 mm

**Ø 4.0**

inserting length (L):

				
8 mm	10 mm	12 mm	14 mm	16 mm

**Ø 4.5**

inserting length (L):

				
8 mm	10 mm	12 mm	14 mm	16 mm


**Ø 5.0**

inserting length (L):

				
8 mm	10 mm	12 mm	14 mm	16 mm

# Conefit implant system

## Sizes available of the Conefit single-phase , one piece „A” type of Implants



The Conefit single-phase implant system consists of implants with five different diameters. Each diameter implant is available in four lengths to be able to find the right solution for each situation. The  $\varnothing 3.5$  diameter implant is exceptionally suitable for thinner bone than average for holding the toothworks in the long run. The diameter of  $\varnothing 3.75$ ,  $\varnothing 4.0$  and  $\varnothing 4.5$  is beneficial for the average bone structure, while the usage of  $\varnothing 5.0$  implants is advantageous for larger bone supply than average.






All our Conefit implants are made of a homogeneous, high-strength titanium alloy.





**Ø 3.5**

inserting length (L):

				
8 mm	10 mm	12 mm	14 mm	16 mm

**Ø 3.75**

inserting length (L):

				
8 mm	10 mm	12 mm	14 mm	16 mm

**Ø 4.0**

inserting length (L):

				
8 mm	10 mm	12 mm	14 mm	16 mm

**Ø 4.5**

inserting length (L):

				
8 mm	10 mm	12 mm	14 mm	16 mm


**Ø 5.0**

inserting length (L):

				
8 mm	10 mm	12 mm	14 mm	16 mm


# Conefit implant system

Sizes available of the Conefit single-phase, one piece „B” type of implants








The Conefit one-piece implant system consists of implants with five different diameters. Each diameter implant is available in four lengths to be able to find the right solution for each situation. The Ø3.5 diameter implant is exceptionally suitable for thinner bone than average for holding the toothworks in the long run. The diameter of Ø3.75, Ø4.0 and Ø4.5 is beneficial for the average bone structure, while the usage of Ø5.0 implants is advantageous for larger bone supply than average.

All our Conefit implants are made of a homogeneous, high-strength titanium alloy.



**Ø 3.5**

inserting length (L):

				
8 mm	10 mm	12 mm	14 mm	16 mm

**Ø 3.75**

inserting length (L):

				
8 mm	10 mm	12 mm	14 mm	16 mm

**Ø 4.0**

inserting length (L):

				
8 mm	10 mm	12 mm	14 mm	16 mm

**Ø 4.5**

inserting length (L):

				
8 mm	10 mm	12 mm	14 mm	16 mm

**Ø 5.0**

inserting length (L):

				
8 mm	10 mm	12 mm	14 mm	16 mm

# The drilling protocol of the **Conefit** implant system

## Drilling Indicators:



- alternatively used



- 1/2 length drilling alternatively




- recommended drilling in 3/4 or in full length

## Ø 3.5 mm drilling protocol of implant



Ø 3.5

implant diameter



colour-coded drill diameter

ø2,0 ø2,5 ø2,8 ø3,2 ø3,7 ø4,2 ø5,2


soft bone  
hard bone

## Ø 3.75 mm drilling protocol of implant



Ø 3.75

implant diameter



colour-coded drill diameter

ø2,0 ø2,5 ø2,8 ø3,2 ø3,7 ø4,2 ø5,2



soft bone  
hard bone

# Ø 4.0 mm

drilling protocol of implant

implant diameter

Ø 4.0





# Ø 4.5 mm

drilling protocol of implant

implant diameter

Ø 4.5




# Ø 5.0 mm

drilling protocol of implant


implant diameter

Ø 5.0



# CONEFIT


## Abutments





## **ABUTMENT SYSTEM | TABLE OF CONTENTS**

Bi-phase abutments of Conefit .....	30
Sizes available of Conefit abutments.....	32
Accessories of abutments .....	36



# The abutments of the **Conefit** bi-phase implants

## 1. Implants

### PROSTHETIC ELEMENTS


2. Technical implant
3. Healing cap, narrow
4. Healing cap, anatomical
5. Sampling head for closed spoon
6. Sampling head for open spoon
7. Temporary head, through-bolted, positioned
8. Temporary head, through-bolted, non-positioned
9. Temporary head, through-bolted, positioned, PEEK
10. Temporary head, through-bolted, non-positioned, PEEK

### FOR GLUABLE TOOTHWORK

11. Narrow head, straight
12. Narrow head oblique
13. Universal head, straight
14. Universal head, oblique
15. Anatomical head, straight
16. Anatomical head, oblique
17. Trapezoidal head
18. Delta head

### IMPLANT-LEVEL CASTING HEADS

19. Castable plastic head, positioned
20. Castable plastic head, non-positioned





- 21. Cobalt chromium-based casting head, positioned
- 22. Cobalt chromium-based casting head, non-positioned

#### FOR REMOVABLE TOOTHWORK

- 23. Ball-head
- 24. Locator head, straight


#### FOR SCREWED TOOTHWORK

- 25. Multi-unit head, straight
- 26. Multi-unit head, through-bolted
- 27. Multi-unit head, oblique
- 28. MC head, oblique
- 29. Multi-unit SR head, screwable
- 30. Multi-unit SR head, through-bolted

- 31. Multi-unit SR head, oblique


#### ELEMENTS OF CAD-CAM SYSTEM

- 32. Titanium base
- 33. Press ceramic base
- 34. Tube-head, positioned
- 35. Tube-head, non-positioned
- 36. Scanbody head, through-bolted
- 37. Scanbody head, screwable
- 38. Scanbody head, through-bolted
- 39. Scanbody head, screwable




# Abutments of the **Conefit bi-phase** implants

TECHNICAL IMPLANT




TEMPORARY HEAD, THROUGH-BOLTED, POSITIONED



HEALING CAP, NARROW




TEMPORARY HEAD, THROUGH-BOLTED, NON-POSITIONED




HEALING CAP, ANATOMICAL




TEMPORARY HEAD, THROUGH-BOLTED, POSITIONED, PEEK




SAMPLING HEAD FOR CLOSED SPOON




TEMPORARY HEAD, THROUGH-BOLTED, NON-POSITIONED, PEEK




SAMPLING HEAD FOR OPEN SPOON




NARROW HEAD, STRAIGHT




NARROW HEAD, OBLIQUE 15°; 25°




TRAPEZOIDAL HEAD 15°; 25°




UNIVERSAL HEAD, STRAIGHT




DELTA HEAD 15°; 25°




UNIVERSAL HEAD, OBLIQUE 15°; 25°; 35°; 45°




CASTABLE PLASTIC HEAD NARROW, NON-POSITIONED




ANATOMICAL HEAD, STRAIGHT




CASTABLE PLASTIC HEAD UNIVERSAL, POSITIONED



ANATOMICAL HEAD, OBLIQUE 15°; 25°



Co-Cr BASED CASTING HEAD, POSITIONED



# Abutments of the **Conefit bi-phase** implants

## Co-Cr BASED CASTING HEAD, NON-POSITIONED

platform (mm):  $\varnothing 3,1$   
 $\odot$   $\odot$   
 platform (mm):  $\varnothing 3,5$   
 $\odot$   $\odot$   $\odot$   
 head height (L mm)  
 12

## MULTI-UNIT HEAD, STRAIGHT

platform (mm):  $\varnothing 3,1$   
 $\odot$   $\odot$   
 platform (mm):  $\varnothing 3,5$   
 $\odot$   $\odot$   $\odot$   
 sulcus height (S mm)  
 0,5 1 2 3 4 5 6 7

## INTERFACE, POSITIONED

platform (mm):  $\varnothing 3,1$   
 $\odot$   $\odot$   
 platform (mm):  $\varnothing 3,5$   
 $\odot$   $\odot$   $\odot$

## MULTI-UNIT HEAD, THROUGH-BOLTED

platform (mm):  $\varnothing 3,1$   
 $\odot$   $\odot$   
 platform (mm):  $\varnothing 3,5$   
 $\odot$   $\odot$   $\odot$   
 sulcus height (S mm)  
 0,5 1 2 3 4 5 6 7

## INTERFACE, NON-POSITIONED

platform (mm):  $\varnothing 3,1$   
 $\odot$   $\odot$   
 platform (mm):  $\varnothing 3,5$   
 $\odot$   $\odot$   $\odot$

## MULTI-UNIT HEAD OBLIQUE, 20°, 30°

platform (mm):  $\varnothing 3,1$   
 $\odot$   $\odot$   
 platform (mm):  $\varnothing 3,5$   
 $\odot$   $\odot$   $\odot$   
 sulcus height (S mm)  
 1 2 3 4 5 6 7

## BALL-HEAD

platform (mm):  $\varnothing 3,1$   
 $\odot$   $\odot$   
 platform (mm):  $\varnothing 3,5$   
 $\odot$   $\odot$   $\odot$   
 sulcus height (S mm)  
 0,5 1 2 3 4 5 6 7

## MC HEAD OBLIQUE, 20°; 30°

platform (mm):  $\varnothing 3,1$   
 $\odot$   $\odot$   
 platform (mm):  $\varnothing 3,5$   
 $\odot$   $\odot$   $\odot$   
 sulcus height (S mm)  
 1 2 3 4 5 6 7


## LOCATOR HEAD, STRAIGHT

platform (mm):  $\varnothing 3,1$   
 $\odot$   $\odot$   
 platform (mm):  $\varnothing 3,5$   
 $\odot$   $\odot$   $\odot$   
 sulcus height (S mm)  
 0,5 1 2 3 4 5 6 7


## MULTI-UNIT SR HEAD, SCREWABLE

platform (mm):  $\varnothing 3,1$   
 $\odot$   $\odot$   
 platform (mm):  $\varnothing 3,5$   
 $\odot$   $\odot$   $\odot$   
 sulcus height (S mm)  
 0,5 1 2 3 4 5 6 7


MULTI-UNIT SR HEAD, THROUGH-BOLTED




TUBE HEAD, POSITIONED




MULTI-UNIT SR HEAD, OBLIQUE 20°; 30°




TUBE HEAD, NON-POSITIONED




TITANIUM BASE




SCANBODY HEAD, THROUGH-BOLTED



PRESS CERAMIC BEASE








SCANBODY HEAD, SCREWABLE



# Accessories of the **Conefit abutments**






## HEADS CREW, SURGICAL



platform (mm):  $\varnothing 3,1$   
   
 platform (mm):  $\varnothing 3,5$   
  






## MULTI-UNIT CUP SCREW



platform (mm):  $\varnothing 3,1$   
   
 platform (mm):  $\varnothing 3,5$   
  






## SAMPLING HEAD SCREW, FOR CLOSED SPOON



platform (mm):  $\varnothing 3,1$   
   
 platform (mm):  $\varnothing 3,5$   
  

## SR-HEADSCREW








platform (mm):  $\varnothing 3,1$   
   
 platform (mm):  $\varnothing 3,5$   
  

It is the same in every diameter






## SAMPLING HEADS CREW, FOR OPEN SPOON



platform (mm):  $\varnothing 3,1$   
   
 platform (mm):  $\varnothing 3,5$   
  






## SR-THROUGH-BOLT



platform (mm):  $\varnothing 3,1$   
   
 platform (mm):  $\varnothing 3,5$   
  

## MULTI-UNIT HEAD SCREW








platform (mm):  $\varnothing 3,1$   
   
 platform (mm):  $\varnothing 3,5$   
  

It is the same in every diameter


## INTERFACE SCREW








platform (mm):  $\varnothing 3,1$   
   
 platform (mm):  $\varnothing 3,5$   
  

It is the same in every diameter





## MULTI-UNIT THROUGH-BOLT



platform (mm):  $\varnothing 3,1$   
   
 platform (mm):  $\varnothing 3,5$   
  





BALL-HEAD CAP. NORMAL



platform (mm):  $\varnothing 3,1$   
  
 platform (mm):  $\varnothing 3,5$   
    
 sphere diameter 2,5 mm  
 It is the same in every diameter





CASTABLE HEAD, PLASTIC



platform (mm):  $\varnothing 3,1$   
  
 platform (mm):  $\varnothing 3,5$   
    
 It is the same in every diameter





BALL-HEAD CAP. MICRO



platform (mm):  $\varnothing 3,1$   
  
 platform (mm):  $\varnothing 3,5$   
    
 sphere diameter 1,8 mm  
 It is the same in every diameter

Co-Cr BASED CASTABLE HEAD



platform (mm):  $\varnothing 3,1$   
  
 platform (mm):  $\varnothing 3,5$   
    
 It is the same in every diameter





LOCATOR HEAD CAP SET



pink cap:  
 10-20° deviations, 3lbs retention  
 blue cap:  
 10-20° deviations, 1,5lbs retention  
 red cap:  
 20-40° deviations, 1lbs retention  
 It is the same in every diameter.





CASTABLE HEAD FOR SR-HEAD



platform (mm):  $\varnothing 3,1$   
  
 platform (mm):  $\varnothing 3,5$   
    
 It is the same in every diameter





INTERFACE PLASTIC




platform (mm):  $\varnothing 3,1$   
  
 platform (mm):  $\varnothing 3,5$   
    
 It is the same in every diameter

Co-Cr BASED CASTABLE HEAD FOR SR-HEAD



platform (mm):  $\varnothing 3,1$   
  
 platform (mm):  $\varnothing 3,5$   
    
 It is the same in every diameter

**CONEFIT**  
Instrument kits








## INSTRUMENT KITS | TABLE OF CONTENTS

Large instrument kit .....	40
Small instrument kit .....	40
Surgical drills .....	44
Ratchet torque wrench .....	46




# Conefit Large instrument kit

Our set of instruments includes instruments essential for dental implantation. The Conefit Large Instrument Kit contains 33 instruments, in a wide range of sizes, for a wide range of applications. The trays are structured


according to the surgical order, labeling makes their usage easier. The tray is also suitable for sterilizing the instruments, this can be done together with the box or separately, as the tray can be uplifted from the box.

The layout of the elements of **Conefit Large instrument kit**



**Preparation**

- 1. Spear-pointed drill
- 2. Pre-drill Ø2.0
- 3. Pre-drill Ø2.5
- 4. Pre-drill Ø 2.7
- 5. Depth gauge




**Thread formation**

- 6. Core drill, Ø 2.8
- 7. Core drill, Ø 3.2
- 8. Core drill, Ø 3.7
- 9. Core drill, Ø 4.2
- 10. Core drill, Ø 4.7




**Thread calibration**

- 11. Thread calibrator drill Ø 3.2
- 12. Thread calibrator drill Ø 3.7
- 13. Thread calibrator drill Ø 4.2
- 14. Thread calibrator drill Ø 4.7
- 15. Thread calibrator drill Ø 5.5




**Paralleling pins**

- 16. Paralleling pin, slim
- 17. Paralleling pin, thick




**Implant drivers, mechanical**

- 18. 6LT 2,1 x L6
- 19. 6LT 2,1 x L12
- 20. 6LT 2,4 x L6
- 21. 6LT 2,4 x L12




**Implant drivers, manual**

- 22. 6LT 2,1 x L6
- 23. 6LT 2,1 x L12
- 24. 6LT 2,4 x L6
- 25. 6LT 2,4 x L12



**Prosthetic tools**

- 26. Gingiva punch Ø4,0xL15
- 27. Head wrench, mechanical 6Lt1,27xL10
- 28. Head wrench, mechanical 6Lt1,27xL15
- 29. Head wrench, manual 6Lt1,27xL10
- 30. Head wrench, manual 6Lt1,27xL15
- 31. Locator key
- 32. Multi-unit key 6LT2,7



**Ratchet torque wrench**


- 33. Ratchet torque wrench



# Conefit Small instrument kit


The Conefit Small Instrument Kit is a practical, more economical solution. The instruments that are inevitable for surgeries can be found in the small tray as they can be found in the large tray as well, only in a smaller range of sizes. There are 19 instruments in the small instrument tray.

The layout of the elements of **Conefit Small instrument kit**




### Preparation

- 1. Spear-pointed drill
- 2. Pre-drill Ø2.0
- 3. Pre-drill Ø2.5
- 4. Depth gauge




### Thread formation

- 5. Core drill, Ø 2.8
- 6. Core drill, Ø 3.2
- 7. Core drill, Ø 3.7
- 8. Core drill, Ø 4.2
- 9. Core drill, Ø 4.7




### Paralleling pins

- 10. Paralleling pin, slim
- 11. Paralleling pin, thick




### Implant drivers

- 12. Mechanical, 6LT 2,4 x L12
- 13. Manual, 6LT 2,4 x L12



### Prosthetic tools

- 14. Gingiva punch Ø4,0xL15
- 15. Head wrench, mechanical 6LT1,27xL10
- 16. Head wrench, manual 6LT1,27xL10
- 17. Locator key
- 18. Multi-unit/ball-head key 6LT2,7




### Ratchet torque wrench

- 19. Ratchet torque wrench



# Conefit surgical drills



**The BIONIKA drills** - which can be used during implant insertion - are available in a wide range of sizes (compatible with different instrument kits to provide you with the most economical solution).

Our drills are externally cooled and have bone collecting properties. Acidic alloy steel and excellent sharpness guarantee long-term use. Each drill is provided with the required drilling depths. Diameters are indicated by color codes.

**Spiral drills** are suitable for preparing implant nests. They are recommended to use according to the drilling protocol, in the case softer and harder bone structures.

**The Thread Calibrator Drills** are suitable for expanding the implant nests as needed, so that we can extend the upper third of the bone nest. They are recommended to use in the case of harder than average bone structure.

## Sizes available of the Conefit surgical drills

### Pre-drills

Ø2.8



Ø3.2



Ø3.7



### Core drills

Ø2.8



Ø3.2



Ø3.7



Ø4.2



Ø5.2



### Thread groove calibrator drills

Ø3.2



Ø3.7



Ø4.2



Ø4.7




Ø5.2



# Ratchet torque wrench

The ratchet torque wrench is used to tighten and insert screws and implants. Using pre-set torque, this prevents the implant from fracturing and ensures the optimum power transfer when inserting the implant. The scale of the torque rates from 15 to 35 Ncm. The desired torques can be adjusted from 15 Ncm to the right for the desired values.

The torque of the key can be infinite if the adjusting screw of ratchet torque wrench screwed to the stop, thus it can also be used for producing a much greater torque than the torque illustrated on the scale if its needed.





















- ERGONOMIC CONSTRUCTION
- RELIABLE LIFETIME
- PRECISE TORQUE CONTROL
- STAINLESS STEEL
- SOUND SIGNAL
- REPLACABLE KEY INTERLINE
- STABLE TORQUE SCALE DURING USE
- PROFESSIONAL QUALITY

- 15 Ncm
- 25 Ncm
- 35 Ncm
- endless torque

The desired torques can be adjusted from 15 Ncm to the right.



# Applications of Ratchet torque wrench

Heads and Screws	Key Interline	Torque	
<b>Locking screw</b>			
Healing cap			
<b>Sampling head screw</b> for closed and open spoon		<p><b>Manual key driver</b> 10-15 Ncm</p>	
<b>Sampling head</b> for closed and open spoon			
<b>Head screw</b>			
Universal head, straight			<p><b>Ratchet torque wrench</b></p> <p>Torque of the required screw tightening:</p> <p>In the case of M1,4 screw it is 15 Ncm</p> <p>In the case of M1,6 screw it is 20 Ncm</p> <p>In the case of M1,8 screw it is 25 Ncm</p> <p>In the case of M2,0 screw it is 30 Ncm</p>
Universal head, oblique			
Anatomical head, straight			
Anatomical head, oblique			
Titanium base			
Multi-unit head, through-bolted			
<b>Multi-unit head screw, SR-head screw</b>			
Multi-unit head, screwable			
SR-head, screwable			
Ball-head			
Locator head			



# Get to know our other products!

Ask for our catalogs or visit one of our websites below:



[www.bionika.hu](http://www.bionika.hu)



[www.shop.bionika.hu](http://www.shop.bionika.hu)

